VH Suction Pump ***Tabletop Aspirator System***

Model 11-100

Model 11-110


**Please read this instruction manual
completely before operating this unit.**

INSTRUCTION MANUAL • ENGLISH & ESPAÑOL

Stop!
Please ensure you have all of the following components before using your device.

Suction Pump
Motor Housing

Air Filter
2x 4" Tube

6' Patient Tubing
(mouthpiece catheter
NOT included)

Collection Jar
with Lid


(5) Air Filters

2-pc Power Cord

Auto Adapter
(11-110 only)

Instruction Manual
Replacement Parts
Order Form

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

INDEX

Care & Safety Information	4-5
Introduction	6
Indications for Use.....	6
Product Features.....	7
Set Up & Operation	8-9
Battery/Power Operation	10
Cleaning and Disinfection	11-12
Changing the Filter	13
Troubleshooting	14
Warnings & Cautions.	14
Device Symbols Warnings	15
Electromagnetic Compatibility	16-19
FCC Statement.....	20
Product Specifications	21
Warranty	22
Instrucciones en Español	23-44

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1103 08/18
©2018 Veridian Healthcare, LLC

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product. Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your suction pump may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the nebulizer's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Never allow the suction pump to be used without adult supervision.
- Do not use the suction pump for anything except patient treatment as directed by a physician. Do not use the suction pump in a manner inconsistent with these instructions.
- The suction pump is an electronic device, keep away from water to avoid electrical shock. Do not immerse the unit in liquid and do not use while bathing. If the unit should fall into water, do not reach for it without first unplugging the unit.
- Do not use while driving.
- Unplug the unit before cleaning it and after each use. Do not leave the unit plugged in when not in use.
- Use only the included power cord(s).
- Always ensure that the air vents are not blocked or obstructed in any way.
- Dropping the suction pump or causing other heavy impact may damage the unit and affect proper operation.
- Do not expose the unit to direct sunlight, humid environments, extreme temperatures, strong static electricity or electromagnetic waves.

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

- If the suction pump is not functioning correctly, if it makes unusual sounds, or if during use the patient experiences pain or discomfort, immediately stop use.
- The patient tubing and catheter should be cleaned daily and replaced monthly during regular use or at the prescribing physician's instructions.
- Do not use the unit if it has been dropped to the ground, exposed to extremetemperatures or high humidity, or damaged in any way.
- Do not disassemble the suction pump or attempt to repair it.

INTRODUCTION




Thank you for purchasing Suction Pump Aspirator. Please read this instruction manual completely before using your unit to ensure safe and effective use, as well as long life of the product.

WARNING: Federal law allows suction pumps to be sold by a physician's prescription only. These devices should be used only at the direction of a doctor or healthcare professional and should be used only as intended and described in this manual.

INDICATIONS FOR USE

The Suction Pump Tabletop Aspirator System is intended to be used to remove infectious materials from wounds or fluids from a patient's airway or respiratory support system.

PRODUCT FEATURES


SET UP & OPERATION

Before the first use of the suction pump if the unit has not been used for an extended period of time, it should be cleaned thoroughly as described in the Cleaning and Disinfection section of this manual.

Assemble the tubing and filter set:

- Firmly push the 4" tubing onto each end of the filter. The tubing should fit tightly onto the filter. If the tubing is not snug, then the set should be replaced.
- Tubing/filter set should be left unassembled when not in use to prevent any moisture build up in the tube or filter.


Assemble the tubing and filter set:

1. Replace the lid on the collection canister. The jar lid will not lock or fasten to the canister but should fit snugly and securely; take care when transporting liquid.
2. Place the canister into the open ring on the suction pump body.

3. Attach the tubing/filter set to the suction pump body, then securely attach the 6' length patient tubing to the port on the collection canister lid.

SET UP & OPERATION

Always inspect the tubing and mouthpiece (not included) for cracks or wear prior to use and replace as necessary.

Patient use:

1. Place the suction unit on a flat, stable surface within operating range of the patient and an electrical outlet. Ensure that the air vents on the side of the unit housing are not blocked or obstructed.
2. Plug the power cord into the electrical outlet.
3. Press the power switch to turn the unit on and activate the suction operation.
4. Verify the suction unit is at the desired suction level prior to patient suction by turning the adjustment knob until the required level is reached. The prescribing physician can advise and explain what level is necessary for the specific patient's care.
NOTE: This device does not include the patient catheter (commonly called suction tip or pipette) as this piece is meant to be replaced frequently. Consult your physician for advice and information on obtaining the patient catheter.
5. When complete, turn the unit off and unplug it from the outlet if it will not be used again for a period of time.
6. Dispose of the contents following the directions of the prescribing physician.
7. Follow the Cleaning and Disinfection directions as advised in this manual.

BATTERY/ POWER OPERATION

The Suction Pump Tabletop Aspirator System is offered with two power option models — 2 piece DC Power Adapter; model 11-110 is rechargeable and includes an auto adapter.

Use common sense electric device caution when operating this unit. Always turn off the power and unplug the device before cleaning or if the unit is not working properly.

- The included power cord should be plugged into the back of the suction pump unit.
- unplug the power cord from the wall and unit when the suction pump is not in use.

Rechargeable DC Power (11-110)

- This unit can operate by the internal rechargeable battery or under AC power.
- The included 2-piece power cord can be plugged into the unit for use or to charge the device. The auto-adapter is intended for charging the internal battery on the go; do not use this product while driving.
- Assemble the power cord and plug it into the unit prior to plugging into a standard wall socket. Operating the unit while plugged into a wall socket will not affect the charged battery life remaining on the internal battery.
- Charge time from expired battery: 8 hours.
- Charged battery life: 60 minutes at normal use; 50 minutes at maximum vacuum flow.

CLEANING & DISINFECTION

The patient tubing and mouthpiece (not included) should be thoroughly cleaned after each use and disinfected each day when being used frequently unless a different cleaning regimen is specified by your healthcare professional; additional hygienic precautions may be required for use with high-risk patients. Always follow the direction of your healthcare provider.

Cleaning After Each Use:

1. Ensure that the power cord is disconnected from the power outlet.
2. Disconnect the air tube, filter, tube set and mouthpiece (not included). Remove the collection canister and dispose of its contents. Do not wash or submerge the suction pump body.
3. Wash and rinse all pieces in hot water, using a mild, non-abrasive detergent.
4. Allow all pieces to air dry completely, or dry them thoroughly with a clean, soft towel. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.
5. Occasionally wipe the body of the unit with a damp cloth and allow it to air dry. Only clean the body of the unit when it is unplugged and completely cool.

Never use soap pads or powder cleanser to clean any part of the unit. Never use any cleaning products that may be harmful when swallowed or inhaled.

It is recommended that air filter and tubing and collection canister be replaced regularly depending on the frequency of therapy treatments or when damaged, cracked or holes appear in the tubing. To order replacement parts, use the included order form, call Customer Care at 866-326-1313 or contact your local dealer.

CLEANING & DISINFECTION

Disinfection at the End of Day / Treatment Period:

The suction pump should be disinfected once daily, after the last use of the day, unless a different disinfection regimen is specified by a healthcare professional. Disinfect the suction pump as follows:

1. Ensure that the power cord is disconnected from the power outlet.
2. Disconnect the air tube, filter, tube set and mouthpiece (not included). Remove the collection canister and dispose of its contents. Do not wash or submerge the suction pump body.
3. Wash all parts in hot water, using a mild, non-abrasive detergent.
4. Submerge all parts in a solution of one part white vinegar and three parts distilled water for a minimum of 30 minutes.
5. Rinse all pieces thoroughly with water.
6. Allow all pieces to air dry completely, or dry them thoroughly with a clean, soft towel. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.

CHANGING THE FILTER

- Never operate the unit without a filter.
- Under regular use, replace the filter monthly, or sooner if it turns gray, appears contaminated is blocked or becomes wet.
- Use only Veridian Healthcare filters; additional filters are available for order using the included order form, calling Customer Care at 866-326-1313 or contacting your local dealer.
- Do not substitute cotton or any other materials for the filters.
- Do not attempt to clean the filters or reuse a filter.

To replace the filter, follow these steps:

1. Change the filter only when the unit is off, unplugged and completely cool.
2. Gently pull the filter cap out of the compressor.
3. Remove the filter from the filter cap.
4. Place a new filter in the filter cap.
5. Return the filter cap to its place, ensuring that it is firmly in place.

TROUBLESHOOTING

The unit will not turn on.

Check all power connections and ensure the cord is firmly plugged into the unit and wall socket. Allow the rechargeable unit sufficient time to recharge or operate with the power cord plugged into the wall socket.

The motor is running, but there is no vacuum pressure.

Check that the tube connections are secure and that there are no leaks or cracks in the tubing. Increase the vacuum level with the pressure adjustment dial.

There is very little vacuum pressure.

Increase the vacuum level with the pressure adjustment dial. Check that the tube connections are secure and that there are no leaks or cracks in the tubing.

**FOR ASSISTANCE IN TROUBLESHOOTING ADDITIONAL ISSUES, PLEASE CONTACT
CUSTOMER CARE AT 866-326-1313**

WARNING AND CAUTIONS

Suction Pump Aspiration is intended for use as prescribed by physician for an individual patient in medical applications only.

Use common sense electrical caution when operating the device. The VH Suction Pump Tabletop Aspirator System runs on AC power or DC rechargeable battery. Always turn off the power and unplug the device before cleaning or if the unit is not working properly.

Never attempt to disassemble the device or service a non-working device. Contact Veridian Healthcare Customer Care department at 866-326-1313 for assistance in warranty repair and service on your unit.

DEVICE SYMBOLS WARNINGS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use.
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
	Keep Dry—This device should be kept dry; never submerge the unit. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor.
	Double Insulation—This unit is double-protected and thus to protection class 2.

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY

Table 2

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

FCC STATEMENT

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesirable operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

CAUTION: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

PRODUCT SPECIFICATIONS

Name	VH Suction Pump Tabletop Aspirator System
Model Number	11-100, 11-110 (DC rechargeable battery)
Sound Level	Maximum 60 dBA at 1 meter
Maximum Negative Pressure	530 mmHg
Negative Pressure Range	150-530 mmHg
Maximum Suction Speed	25 l/min
Collection Canister Capacity	800 cc
Operating Conditions	Temperature: 50°F – 104°F (10°C – 40°C); Humidity: 0% – 95%
Storage Environment	Temperature: -13°F – 158°F (-25°C – 70°C); Humidity: 0% – 95%
Unit Dimensions (W x H x D)	14.3" x 8.4" x 6.9"
Unit Weight	3.75 lbs (11-100); 4.52 lbs (11-110)
Power Connection	12 V
Power Consumption	48 VA Max
Included Accessories	Suction pump unit, plastic collection canister with lid, 6' patient tubing, 2x 4" tubing/filter set, five air filters, 2 -piece DC Power Adapter, Auto Adapter (11-110)

Specifications are subject to change without notice.

2-YEAR LIMITED WARRANTY

Veridian Healthcare guarantees that its VH Suction Pump Tabletop Aspirator System will be free from manufacturing defects under normal use for a period of two-years. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty is applicable to the suction pump unit and motor and does not cover batteries or other power sources that may be provided with or used with the suction pump unit or the collection canister. This warranty is voided if the suction pump unit is misused or abused in any manner.

If the suction pump unit fails to operate within the warranty period while under the original purchaser owns the product, return it postage prepaid to: **Attn: Repair Department, 1175 Lakeside, Gurnee, IL, 60031**; include your name, address, phone number and brief description of the problem with your nebulizer. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. Veridian Healthcare will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you. This warranty constitutes the warrantor's only reasonability and obligation to repair and/or replace materials or components, or refund the purchase price. Veridian will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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