
Product Overview

How to Set Your Pulse Ring

Hold down the SET button for 2 seconds to start setting segments.

1. ON/OFF the clock function
 - Press SET button to select
 - Press MODE button to next
 If "OFF" is selected, move to step 3
 2. Minute and hour
 - Press SET button to advance numbers
 - Press MODE button to next
- Note: Time displayed in military time.

3. Count-up time limit
 - Press SET button to select time limit
 - Press MODE button to next
 Ex. 00:30 indicates that the stopwatch will automatically stop after 30 minutes from starting.
4. Minimum target heart rate limit
 - Press SET button to set your minimum heart rate limit
 - Press MODE button to next
5. Maximum target heart rate limit
 - Press SET button to set your maximum heart rate limit
 - Press MODE button to next

Setting is now complete. Press MODE button to exit settings.


How to Monitor Your Heart Rate

1. Press MODE button to turn on measurements.
 2. The monitor starts to catch the signal and "-- --" and MODE flash.
 3. The current heart rate with a flashing will display.
- To exit, press the MODE button again.


While the heart rate measurement is below or over the target, "MIN" or "MAX" flashes with alarm beeps.

For the best performance, wear your Pulse Ring on the finger with the best fit. Keep it from moving or loosening while taking measurements.

You may attach the enclosed finger pad onto the battery cover for optimal fit.

Error Message


"Err" flashes with alarm beeps when the signal is not caught within 30 seconds. "OFF" will display on the screen, indicating the heart rate monitoring function is stopped. Please refit your Pulse Ring and start again.

How to Use the Stopwatch

1. Press SET button to activate stopwatch.
 2. Press the SET button to start counting.
 3. Press the SET button to stop counting.
- To resume, press MODE button once. To exit, press MODE button twice.

Battery Replacement


One 3V CR2032 lithium battery is needed. Replace the battery when the color of reading display becomes lighter, or the monitor has no reaction by pressing any button.

Maintenance

We suggest wiping your Pulse Ring with furniture polish (wax) to refresh the unit.

Using the Case and Lanyard

Target Heart Rate Zone

To find your maximum heart rate, you can use the age formula, max HR = 220 - age.

Maximum heart rate, or max HR, is the highest number of heartbeats per minute (bpm) in an all-out effort. Training intensities can be expressed as percentages or max HR. There are three different exercise zones: hard, moderate and light.

HARD

Benefits: maximizes performance capacity
 Recommended for: fit people and short exercise sessions

MODERATE

Benefits: improves aerobic fitness
 Recommended for: everybody and for sessions of moderate length

LIGHT

Benefits: manages weight and general wellness
 Recommended for: everybody, for longer sessions

AGE	maxHR 100%	LIGHT 60-70%	MODERATE 70-80%	HARD 80-90%
20	200	120 - 140	140 - 160	160 - 180
30	190	114 - 133	133 - 152	152 - 171
40	180	108 - 126	126 - 144	144 - 162
50	170	102 - 119	119 - 136	136 - 153
60	160	96 - 112	112 - 128	128 - 144
beats per minute				

IMPORTANT! Please be aware that these are general guidelines. Target zones will vary from person to person. For goal settings, please consult a trained professional.

Specifications

- Heart rate monitor 30 - 250 beats per minute
- Stopwatch 00:00:00 to 10:00:00 (hrs:min:sec)
- Real time clock 00:00:00 to 23:59:59 (hrs:min:sec)
- Operating Temperature 32°F - 122°F (0°C - 50°C)
- Storage Temperature -4°F - 140°F (-20°C - 60°C)
- Battery One 3V CR2032 lithium battery
- Contents Pulse Ring (battery inside), case, lanyard, extra battery, instructions

Your heart rate ring is covered by a 30-day warranty. For details of the warranty please contact Customer Care at 800-622-4714.

This is not a medical device and can not be used for medical diagnostic and treatment activities.

Manufactured for:
 Briggs Medical Service Company
 1931 Norman Drive
 Waukegan, IL 60085

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©2011 Briggs Medical Service Company