
20-Second Flexible Tip Dual Scale Digital Thermometer

Thank you for purchasing a Digital Thermometer. Please read the following detailed instructions to achieve the most accurate temperatures and safe operation. With proper care and use, your thermometer will provide you with years of reliable readings.

What is a “Normal” Temperature?

Although the generally accepted “normal” temperature reading is 98.6°F (37.0°C), temperature readings can vary from 97°F (36.1°C) to 99°F (37.2°C) and still may be considered “normal”. Variations in temperature can be attributed to activities such as exercise, smoking, eating and drinking. Even time of day may influence your temperature. For example, your temperature is lower in the morning than in the afternoon. Other variations may be due to the location of the temperature reading. While the oral temperature follows the guidelines previously described, a rectal temperature is generally 1.0°F (0.5°C) higher. Conversely, an axillary (under the arm) temperature will be 1.0°F (0.5°C) lower.

Features of your MABIS Digital Thermometer BATTERY COVER (Animal May Vary)

Cleaning Instructions

You can clean your Digital Thermometer by wiping it with soap and water or 70% Isopropyl Alcohol. THE UNIT MUST NEVER BE SUBMERSED IN ANY LIQUID OR STERILIZED IN ANY WAY (GAS OR STEAM AUTOCLAVES). Failure to follow these instructions will cause severe damage to your unit and void your warranty.

Switching Dual Scale

NOTE: To change between Celcius and Fahrenheit readings, hold the POWER button for 2 seconds until the “Lo” with an “m” displays and changes between °F and °C. Press the POWER button again to turn the thermometer off.

How to Take Your Temperature

- 1. Clean the thermometer according the instructions provided.
- 2. Apply a disposable probe cover to the clean, dry thermometer.
- 3. Determine whether the measurement site will be oral, rectal or underarm.

NOTE: Measurement time may vary based on the selected site.

100.0°F

Lo

98.6°F

Lo

Once the desired temperature is selected, press the POWER button (the unit will beep) to activate the thermometer. This is the LCD display test, which insures all parts of the display are functioning properly.

The last measured temperature with an “m” will appear. If there is no previous measurement, “Lo” °F(°C) will be displayed.

98.6°F(37.0°C) will briefly appear

Finally, a “Lo” with a flashing °F(°C) will appear in the display to indicate that the thermometer is ready for measurement.

NOTE: If the room temperature is above 89.6°F(32.0°C), the room temperature will display instead of the Lo.

Taking an Oral Temperature Approximate measurement time: 20 seconds

1. Using the diagram on the right, place the probe tip under the tongue as near as possible to the heat pocket. The diagram also shows the temperature variables by not placing the tip in the heat pocket area.

NOTE: Drinking hot or cold fluids, exercising, smoking or performing other activities may raise or lower your temperature. Therefore, it is important to relax for approximately 5 minutes with your mouth closed prior to taking a reading.

2. The current temperature is shown and the °F (°C) will continue to flash during the reading.

3. When the peak temperature has been reached, the °F(°C) will stop flashing and the thermometer will sound ten long beeps if the temperature is below 100°F(37.8°C). If the temperature is above 100°F(37.8°C) the unit will sound ten sets of three short beeps each.

4. To turn the unit off, press the ON/OFF button. The unit will automatically shut off after approximately nine minutes.

NOTE: The mouth must remain closed during the measurement for the most accurate reading. Opening the mouth could result in an extended measurement time and may affect the reading.

Taking an Axillary Temperature (under the arm) Approximate measurement time: 30 seconds

- 1. Wipe the underarm with a dry towel.
 - 2. Place the probe tip under the arm so the tip is touching the skin with the thermometer perpendicular to the body.
 - 3. Position the arm across the chest so the probe tip is well covered by the arm. This also insures the probe is not affected by the room's air.
 - 4. When the peak temperature has been reached, the °F(°C) will stop flashing and the thermometer will sound ten long beeps if the temperature is below 100°F(37.8°C). If the temperature is above 100°F(37.8°C) the unit will sound ten sets of three short beeps each.
 - 5. To turn the unit off, press the ON/OFF button. The unit will automatically shut off after approximately nine minutes.
- NOTE: Axillary temperatures are generally 1.0°F (0.5°C) lower than oral readings.

Taking a Rectal Temperature Approximate measurement time: 15 seconds

- NOTE: Rectal temperatures are an appropriate and reliable method for infants and small children. Otherwise, this method should only be used when it is impossible or impractical to take an oral or axillary temperature.
- 1. After applying the disposable probe cover, apply a water-soluble lubricant to the tip.
 - 2. Gently insert the probe a MAXIMUM of 1/2” into the rectum.
 - 3. When the peak temperature has been reached, the °F(°C) will stop flashing and the thermometer will sound ten long beeps if the temperature is below 100°F(37.8°C). If the temperature is above 100°F(37.8°C) the unit will sound ten sets of three short beeps each.
 - 4. To turn the unit off, press the ON/OFF button. The unit will automatically shut off after approximately nine minutes.
- NOTE: Rectal temperatures are generally 1.0°F (0.5°C) higher than oral readings.

Using Disposable Probe Covers

It is recommended to use the enclosed disposable probe covers to prevent cross contamination. Additional probe covers can be purchased from your Dealer or with the order form provided.

General Care

- 1. The thermometer should only be used under the supervision of an adult.
 - 2. Do not walk, run or talk during temperature taking.
 - 3. Clean the thermometer before and after each use.
 - 4. Store the unit in the protective case when not in use.
 - 5. Do not bite the probe.
 - 6. Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
 - 7. Dropping or subjecting your thermometer to strong shocks should be avoided.
 - 8. Do not attempt to disassemble the unit, except to replace the battery.
 - 9. If measurement is not attainable, the unit must be replaced.
- NOTE: Performance of the device may be degraded if: operated outside of stated temperature and humidity range; stored outside of stated temperature and humidity range; thermometer undergoes mechanical shock (drop); patient temperature is below ambient temperature.

Changing the Battery

The battery in your digital thermometer needs to be replaced when the “ ” symbol appears. Replace the battery with a LR41 or equivalent.

To replace the battery,

1. Remove screw and slide the battery cap away from the unit.
2. Gently pull the battery tray out until it stops. Using a non-metal pointed tool, carefully remove the battery and insert the new battery with the positive (+) side facing up.
3. Carefully slide the battery cap back into place, making sure that the cap fits snugly to protect the battery compartment from moisture and replace screw.

NOTE: Please properly dispose of the battery away from children and heat.

Specifications

DISPLAY:	Liquid Crystal Display	
DISPLAY RESOLUTION:	4 Digits (1/10TH of a degree increments)	
MEASUREMENT RANGE:	90.0°F to 109.9°F (32°C - 42.9°C)	
ACCURACY:	±0.2°F (98.0°F - 102.0°F) ±0.1°C (37.0°C - 39.0°C)	
AVERAGE MEASUREMENT TIME:	Oral 20 seconds. Rectal 15 seconds. Axillary 30 seconds	
tone:	Thermometer beeps when activated and deactivated. A series of beeps signal when peak temperature is reached.	
MEMORY:	Displays last measured temperature	
BATTERY:	One, LR41 or equivalent	
BATTERY LIFE:	Approximately 1,500 measurements	
AUTOMATIC SHUT-OFF:	After approximately 9 minutes	
AMBIENT OPERATING TEMPERATURE:	50.0°F to 113.0°F (10°C to 45°C)	
HUMIDITY:	15% to 95% RH	
STORAGE TEMPERATURE:	-4.0°F to 122.0°F (-20°C to 50°C)	
HUMIDITY:	15% to 95% RH	
DIMENSIONS:	Puppy: 5-3/4" x 7/8" x 3/4" Bunny: 5-7/8" x 7/8" x 5/8" Duck: 5-3/4" x 7/8" x 3/4"	Frog: 5-3/4" x 1-1/8" x 3/4" Cow: 5-5/8" x 7/8" x 3/4"
WEIGHT:	Puppy: 0.48 ounces (13.5 grams) Bunny: 0.42 ounces (11.8 grams) Duck: 0.50 ounces (14.2 grams)	Frog: 0.51 ounces (14.6 grams) Cow: 0.47 ounces (13.2 grams)

LIMITED FIVE YEAR WARRANTY

The warrantor guarantees that its Digital Thermometer will be free from manufacturing defects under normal use for a period of five years.

This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Thermometer. This warranty is voided if the Digital Thermometer product is misused or abused in any manner.

If the Digital Thermometer fails to operate during the time the original purchaser owns it, return it postage prepaid to: Repair Department, 1931 Norman Drive South, Waukegan, IL 60085. The warrantor will repair or replace the defective unit at the option of the warrantor. Repair or replacement of the defective unit, at warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantors only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

CUT OUT AND MAIL TO:

WARRANTY REGISTRATION DEPARTMENT
1931 NORMAN DRIVE SOUTH
WAUKEGAN, IL 60085

THERMOMETER WARRANTY REGISTRATION
(MUST BE COMPLETED WITHIN 10 DAYS OF PURCHASE)

NAME: _____

ADDRESS: _____

CITY: _____ STATE: _____ ZIP: _____

PURCHASED FROM: _____

CITY: _____ STATE: _____ ZIP: _____

DATE PURCHASED: _____

MODEL PURCHASED: _____

USER INFORMATION (OPTIONAL)
☐ MALE ☐ FEMALE AGE: _____

AGE OF YOUNGEST FAMILY MEMBER
☐ UNDER 6 MONTHS ☐ 6 MONTHS - 1 YEAR
☐ 1 - 3 YEARS ☐ 3 - 5 YEARS ☐ OLDER

PROBE COVER ORDER FORM

CUT OUT AND MAIL TO:

MABIS HEALTHCARE
1931 NORMAN DRIVE SOUTH
WAUKEGAN, IL 60085

Please send _____ packages of disposable probe covers #15-618-000 (100 probes per package). Enclosed is \$5.00 per package plus \$3.00 for shipping and handling for each order.

PAYMENT METHOD (do not send cash): ☐ Check ☐ Money Order

Please ship to:

NAME: _____

ADDRESS: _____

CITY: _____ STATE: _____ ZIP: _____

PHONE: (_____) _____

NOTE: This order form is for shipments within the United States and Canada.

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